

# MIGHTYmo™ miniMo™

**Digitally powered to allow  
you to confidently and quickly  
speak your mind.**



## Crisp, Clear Communication

Augmented communicators of all physical and cognitive abilities who prefer or require digitized speech can create messages of any length or complexity. The ability to record more than 100 minutes of speech and sound recordings – in addition to those provided with your MightyMo or MiniMo at time of purchase – ensures ample room for needed vocabulary.

## Easy to Use and Setup

Prerecorded page sets provide professionally recorded vocabulary that can be used by children as a starting point for normal daily activities. Creating customized pages with greater meaning for the augmented communicator is easy with MightyMo and MiniMo. Import pictures and sounds to create enhanced visual representations.

## Peace of Mind

Developed specifically for augmented communicators, MightyMo and MiniMo can absorb the impact of daily living. You can rest easy knowing that DynaCare® covers your MightyMo or MiniMo for one full year from the date we ship the product to you.

## Time for the Important Things

By combining the natural message formation strategies for augmented communicators with tools and resources for care teams that reduce the amount of time needed to learn and program the devices, MightyMo and MiniMo allow for greater time spent making meaningful connections with one another and the world.

## SPECIAL FEATURES

### Professionally Recorded Page Sets

Choose from an adult male voice, an adult female voice or a gender-neutral child's voice.

- Preschool Pages
- Adult Pages
- Medical Pages
- Virtual Pages
- Behavior Modification Pages
- Gateway Jumpstart™

### Onboard Universal Remote Control

- Control common household appliances such as TVs, DVD and CD players, and IR phones, etc. with your MightyMo or MiniMo.

### A menu structure and software interface similar to that of the DyaVox DV4 and MT4

- Makes the transition from MightyMo or MiniMo to the more fully featured Series 4 products easy for everyone.

Integrated microphone and auto-leveling help to ensure that messages are recorded and delivered with a high degree of clarity.

Backup to and restore data from your Windows® PC using USB storage device (included with your purchase).



Technology for human expression.

2100 Wharton Street • Suite 400 • Pittsburgh, PA 15203 • 1.866.DYNAVOX • [www.dynavoxtech.com](http://www.dynavoxtech.com)

## ADDITIONAL FEATURES

### Choice of Symbols

- Picture Communication Symbols® (PCS)®
- DynaSyms® (New DynaSyms and Classic DynaSyms)
- Both PCS and DynaSyms are available at an additional charge.

### DynaCare Certified Unlimited Guarantee™

- One-year DynaCare guarantee covers repair or replacement in all instances except theft or loss
- One- to four-year DynaCare guarantee extension available
- For more information, please see our DynaCare Certified unlimited Guarantee information sheet.

## MIGHTYMO AND MINIMO SPECIFICATIONS

### The following items are included with your MightyMo or MiniMo purchase:

One MightyMo or MiniMo speech-output device, power transformer and cord, shoulder strap, five pre-programmed page sets, one symbol set (DynaSyms or PCS), User's Guide, backup CD and USB storage device. Additional accessories are available.

#### DYNAVOX MIGHTYMO SPECIFICATIONS

- **Weight:** 5 lbs 13 oz / 2.6 kg
- **Size:** 12" x 9" x 3" / 30 x 25 x 8 cm
- **Display:** 12.1" / 30.7 cm active matrix color TFT-LCD
- **Battery Information:** 12/3.7 V li-ion cells, pack is 7.4 V @ 13.2 Ahrs
- **Average Battery Run Time:** Continuous use – 12 hours
- **Power Transformer (Charger):** AC 100-240 V  
1.0 A 50/60 Hz DC 12 V 3A
- **Operating System:** WinCE®

#### DYNAVOX MINIMO SPECIFICATIONS

- **Weight:** 3 lbs 5 oz / 1.5 kg
- **Size:** 8.25" x 6" x 2.25" / 21 x 15 x 6 cm
- **Display:** 6.4" / 16.3 cm active matrix color TFT-LCD
- **Battery Information:** 8/3.7 V li-ion cells, pack is 7.4 @ 8.8 Ahrs
- **Average Battery Run Time:** Continuous use – 8 hours
- **Power Transformer (Charger):** AC 100-240 V  
1.0 A 50/60 Hz DC 12 V 3A
- **Operating System:** WinCE®

#### ACCESS AND SELECTION METHODS

- Touch screen with adjustable hold/release times, touch enter/exit modes, visual and audio scanning with auto, switch, two-switch and inverse auto-scan options
- Mouse pause and mouse click
- AutoZoom is available for touch enter selectors

#### SPEECH OUTPUT

- Digitized speech and sound

#### INTERFACES

- Standard switch port – two jacks for standard switches
- 2 USB type B function connectors
- Speaker port for external speaker or private audio scanning
- Integrated microphone and external microphone jack for digitized speech and sound
- Infrared for learning universal remote control

#### OPERATING AND STORAGE ENVIRONMENTS:

- Operating temperature:  
32° to 122° F / 0° to 50° C
- Storage temperature:  
32° to 122° F / 0° to 50° C

#### MODEL NUMBERS

• MightyMo	800047	(Male Voice)
	800046	(Female Voice)
	800044	(Child's Voice)
• MiniMo	800060	(Male Voice)
	800059	(Female Voice)
	800057	(Child's Voice)

#### POPULAR ACCESSORIES

- Keyguards
- PCS or DynaSyms
- Protective Carrying Case
  - User Accessible Carrying Case
- DynaCare Guarantee Extension
- Daessy® compatible mounting plate
- Additional power transformer and cord

**To learn more about these accessories, please refer to the DynaVox Technologies Accessory and Mounting Guides. Or visit us at [dynavoxtech.com](http://dynavoxtech.com)**



**DynaVox**  
TECHNOLOGIES

Technology for human expression.

2100 Wharton Street • Suite 400  
Pittsburgh, PA 15203  
1.866.DYNAVOX (396.2869)  
[www.dynavoxtech.com](http://www.dynavoxtech.com)